A14

STANDING TRANSFER STAND & TURN FRAME – WHEELED WITH BELT

- 1
- Person to shuffle bottom to front of seat.
- Person to lean to one side. Carer to fit belt (elastic section at top) tucking under persons bottom, repeat for other side.
- Secure belt with Velcro straps.
- Carer to extend locking strap fully, and fasten.
- 2
- Move frame in front of person
- Person to position feet on base plate.
- Carer to apply brakes, and hold onto frame providing counterbalance.
- Person to hold frame, lean forward with head raised and stand.
- Once person is standing, thread locking strap through central bar.
- Carer to release brake.
- Ensure person's balance is stable before manoeuvring frame.
- 3
- Person to be positioned so that back of knees touch seat / bed.
- Carer to apply brakes to frame.
- Person to come forward whilst locking strap is released from frame.
- Carer to hold onto frame providing counterbalance whilst person sits holding onto the frame.
- Once seated, carer to remove frame and belt.









