

Community Equipment Assessor & Technician Prescribed Equipment Guide



Bathing Equipment

Bathlift Fixed (Orca) (M68476)	Page 2
Bathlift (Neptune) (H84588)	Page 3
Bathlift Reclining (AquaJoy) (L35749)	Page 4
Bathing Cushion (M18217)	Page 5
Bathlift Replacement Covers (Orca) (M51683)	Page 6
6" Bath Seat (M78419)	Page 7
8" Bath Seat (M78421)	Page 8
12" Bath Seat (M78433)	Page 9
26" Bath Board - Slatted (M78342)	Page 10
27" Bath Board - Slatted (M78354)	Page 11
28" Bath Board - Slatted (M78366)	Page 12
26" Shower Board - Six Slatted (M78389)	Page 13
27" Shower Board - Six Slatted (M78391)	Page 14
28" Shower Board - Six Slatted (M78408)	Page 15
Handle for Shower & Bath Boards (M78378)	Page 16
Shower Stool with Handles (F20524)	Page 17
Shower Chair - Height Adjustable (M48295)	Page 18
Swivel Bath Seat (N65743)	Page 19
Bath Step - Adjustable 2-4" (M11077)	Page 20
Bath Step - 4" Fixed Height (M99515)	Page 21
Portable Shower Screen (N65731)	Page 22

Other Equipment

2" Raised Toilet Seat (N68185)	Page 23
4" Raised Toilet Seat (N68197)	Page 24
6" Raised Toilet Seat (N05266)	Page 25
Toilet Frame (Standard) (M00869)	Page 26
Toilet Frame (Adjustable Width) (M00870)	Page 27
Toilet Frame Floor Fixing Kit (M24308)	Page 28
Toilet Frame with Raised Seat (M66625)	Page 29
Commode (Adjustable Height) (F18864)	Page 30
Trolley (H87668)	Page 31
Buckingham Caddy (L29877)	Page 32
Perching Stool (No Arms & Back) (H87679)	Page 33
Perching Stool (Arms & Back) (F19893)	Page 34
Perching Stool (Arms & Padded Back) (L33236)	Page 35
High Density Cushion (N70854)	Page 36
Leg Lifter with Loop (F19571)	Page 37
Bed Lever (Slatted & Divan) (M48192)	Page 38
Bed Lever Double End (Divan Only) (M54234) ..	Page 39

Information for Prescribers

This document is for use by professionals only and we have indicated next to items, who is able request and order that specific item.

Please do not request provision or the installation* of any of these items from the County Technician Service. All items listed here are available from NRS and should be ordered from them using iRIS.

**Toilet frame fixing kits should be ordered from NRS and a CTS57 raised for a County Technician to fix down.*

Key to Prescribers

TECHNICIAN

This item can be prescribed by KCC's Technician's.

CEAx

This item can be prescribed by Community Equipment Assessors.

Clinical / Practical Reasoning Required

This item needs authorisation and will require clinical / practical reasoning to be given on assessment form.

Last Updated — 10th May 2018

M68476 - Bathlift Fixed Backrest (Orca)

Description

- Side flaps to ease transfers
- FIXED Backrest (does not recline)
- Foldaway Design
- Large seat Depth
- Low Seat Height
- Easy to use hand control
- Two Part Construction

Dimensions

- Seat height at lowest position: 60mm (2½")
- Seat height at highest position: 420mm (16½")
- Seat width with side flaps folded out: 715mm (28¼")
- Seat width with side flaps folded in: 420mm (16½")
- Width of backrest: 360mm (14¼")
- Seat Length: 590mm (23¼")
- Total weight without hand control: 12.1kg (26.7lb)
- Charger weight: 100g (3.5oz)



TECHNICIAN

Maximum user weight: 135kg (21st)

H84588 - Bathlift (Neptune)

Description

- Rechargeable battery powered bathlift
- Secures to bath with 4 x large suckered feet
- Battery monitoring system
- Waterproof hand control
- Position when fitted allows maximum leg room
- Does not recline

Dimensions

- Seat Height: 70mm – 435mm (2¾" - 17")
- Overall Height: 735mm (29")
- Seat Width: 350mm (13¾")
- Seat Width Inc Flaps: 730mm (28¾")
- Seat Depth: 410mm (16")
- Overall Length: 725mm (28½")
- Product Weight: 9.5kg (21lbs)



TECHNICIAN

Clinical / Practical
Reasoning Required

Maximum user weight: 160kg (25st)

L35749 - Bath Lift Reclining Backrest (AquaJoy)

Description

- Reclining back rest
- Waterproof floating hand control
- Rechargeable battery
- Waterproof hand control
- Will not lower if insufficient battery power to raise after bathing

Dimensions

- Seat height lowest: 90mm (3½")
- Seat height highest: 465mm (18¼")
- Height of backrest: 604mm (23¾")
- Seat width: 390mm (15¼")
- Seat width with flaps open: 686mm (27")
- Overall length: 580mm (22¾")
- Weight without battery: 13.6kg (30lbs)



TECHNICIAN

Clinical / Practical
Reasoning Required

Maximum user weight: 170kg (26st)


M18217 - Bathing Cushion (Manger)

Description

- Enables user to be lowered to the bottom of the bath and raises them back out when required.
- User needs reasonable upper body stability.
- An Airflo compressor is supplied to inflate the cushion.
- Waterproof hand control.

Dimensions

- Maximum Height: 400mm (15½")
- Minimum Height: 20mm (¾")
- Seat Width: 580mm (23")
- Seat Depth: 440mm (17½")
- Backrest Width: 540mm (21")
- Backrest Depth 330mm (13")
- Air hose length: 1.8m
- Total product weight: 2kg (4½lbs)



Maximum user weight: 150kg (24st)

TECHNICIAN

Clinical / Practical Reasoning Required

M51683 - Bathlift Cover Set White

Description

- Replacement covers for Orca Bathlift
- Also available in Blue



Covers Only

M55317 - Bathlift Cover Set Blue

Description

- Replacement covers for Orca Bathlift
- Also available in White



Covers Only

TECHNICIAN

M78419 - Bath Seat Slatted 6"

Description

- Slatted bath seat
- Fixes securely to the bath with large suction pads
- Reinforced PVC
- Large slatted seat area for easy drainage
- Can be used with a bath board.

Dimensions

- Height: 150mm (6")
- Width: 457mm (18")
- Depth: 280mm (11")
- Product Weight: 2.7kg (6 lbs)



TECHNICIAN

CEAx

Maximum user weight: 160kg (25st)

M78421 - Bath Seat Slatted 8"

Description

- Slatted bath seat
- Fixes securely to the bath with large suction pads
- Reinforced PVC
- Large slatted seat area for easy drainage
- Can be used with a bath board.

Dimensions

- Height: 200mm (8")
- Width: 457mm (18")
- Depth: 280mm (11")
- Product Weight: 2.8kg (6 lbs)



TECHNICIAN

CEAx

Maximum user weight: 160kg (25st)

M78433 - Bath Seat Slatted 12"

Description

- Slatted bath seat
- Fixes securely to the bath with large suction pads
- Reinforced PVC
- Large slatted seat area for easy drainage
- Can be used with a bath board.

Dimensions

- Height: 305mm (12")
- Width: 457mm (18")
- Depth: 280mm (11")
- Product Weight: 2.9kg (6½lbs)



TECHNICIAN

CEAx

Maximum user weight: 160kg (25st)

M78342 - Bath Board Slatted 26"

Description

- Reinforced PVC slatted bath board.
- Stainless steel fittings.
- Adjustable brackets for secure fitting.
- Optional grab handle to be ordered separately.

Dimensions

- Length: 660mm (26")
- Width: 230mm (9")
- Bracket to end of board: 25 – 90mm (1 – 3½")



TECHNICIAN

CEAx

Maximum User Weight: 190kg (30st)

M78354 - Bath Board Slatted 27"

Description

- Reinforced PVC slatted bath board.
- Stainless steel fittings.
- Adjustable brackets for secure fitting.
- Optional grab handle to be ordered separately.

Dimensions

- Length: 686mm (27")
- Width: 230mm (9")
- Bracket to end of board: 25 – 90mm (1 – 3½")



TECHNICIAN

CEAx

Maximum User Weight: 190kg (30st)

M78366 - Bath Board Slatted 28"

Description

- Reinforced PVC slatted bath board.
- Stainless steel fittings.
- Adjustable brackets for secure fitting.
- Optional grab handle to be ordered separately.

Dimensions

- Length: 711mm (28")
- Width: 230mm (9")
- Bracket to end of board: 25 – 90mm (1 – 3½")



TECHNICIAN

CEAx

Maximum User Weight: 190kg (30st)

M78389 - Shower Board 6-Slatted 26"

Description

- Reinforced PVC slatted shower board.
- Stainless steel fittings.
- Adjustable brackets for secure fitting.
- Optional grab handle to be ordered separately.



Dimensions

- Length: 660mm (26")
- Width: 381mm (15")
- Bracket to end of board: 25 – 90mm (1 – 3½")

TECHNICIAN

CEAx

Maximum User Weight: 190kg (30st)

M78391 - Shower Board 6-Slatted 27"

Description

- Reinforced PVC slatted shower board.
- Stainless steel fittings.
- Adjustable brackets for secure fitting.
- Optional grab handle to be ordered separately.



Dimensions

- Length: 686mm (27")
- Width: 381mm (15")
- Bracket to end of board: 25 – 90mm (1 – 3½")

TECHNICIAN

CEAx

Maximum User Weight: 190kg (30st)

M78408 - Shower Board 6-Slatted 28"

Description

- Reinforced PVC slatted shower board.
- Stainless steel fittings.
- Adjustable brackets for secure fitting.
- Optional grab handle to be ordered separately.



Dimensions

- Length: 711mm (28")
- Width: 381mm (15")
- Bracket to end of board: 25 – 90mm (1 – 3½")

TECHNICIAN

CEAx

Maximum User Weight: 190kg (30st)

M78378 - Handle for Shower & Bath Boards

Description

- For use with slatted bath and shower boards
- Offers additional support when transferring
- Attaches securely to the board
- Can be useful if it is not possible to fit a rail to the wall.

Dimensions

- Length: 190mm (7½")
- Width: 30mm (1¼")
- Height above board: 100mm (4")
- Handhold (dia of bar): 30mm (1¼")
- Handhold width (Inside): 125mm (5")
- Weight 120g (4½ oz)



TECHNICIAN

CEAx

F20524 - Shower Stool with Handles

Description

- Adjustable height shower stool
- Does not rust
- Removable seat
- Does not rust

Dimensions

- Height: 480 - 630mm (19 – 25")
- Width between arms: 485mm (19")
- Seat Size: 380 – 280mm (15 x 11")
- Base size 455 x 380mm (18 x 15")



Please state Height when prescribing this item.

TECHNICIAN

CEAx



Maximum User Weight: 160kg (25st)

M48295 - Shower Chair Height Adjustable

Description

- Shower Chair.
- Height adjustable.
- Integral handles for extra support.
- Drainage holes

Dimensions

- Seat height: 410 – 537mm (16¼ - 21¼")
- Width between arms: 485mm (19")
- Footprint: 545 x 450mm (21 ½ x 17¾")



Please state Height when prescribing this item.

TECHNICIAN

CEAx



Maximum User Weight: 160kg (25st)

N65743 - Swivel Bath Seat

Description

- Swivel bather with arm rests
- Non-slip non-marking rubber ends on frame.
- Adjustable width

Dimensions

- Frame Width external: 660 – 710mm (26 – 28")
- Frame width internal: 500 – 610mm (20 – 24")
- Depth: 550mm (21½")
- Seat Height from top of bath: 130mm (5")
- Armrest height from top of bath: 230mm (9cm)
- Seat width: 460mm (18")
- Seat depth: 420mm (16½")



TECHNICIAN

Maximum User Weight: 160kg (25st)

M11077 - Bath Steps (Adjustable)

Description

- Each step section has non-slip feet so the bath step remains secure.
- The top surface of each section has a textured surface so it's comfortable on the feet.
- Provides excellent resistance to slipping, even when the client has wet feet.
- Support pillars provide extra strength and support to the bath step
- A recess on the side of each section allows the step to be easily picked up
- Manufactured from polypropylene
- Can be cleaned with a mild disinfectant.

Dimensions

- Length: 500mm (19½")
- Width: 395mm (15½")
- Height Options: 51mm (2"), 77mm (3"), 103mm (4")



*****IMPORTANT: CONTAINS
SMALL PARTS, MAY NOT BE
SUITABLE FOR CHILDREN*****

TECHNICIAN

Maximum User Weight: 190kg (30st)

M99515 - Bath Step

Description

- 100mm (4") step for use in bathroom
- Slip-resistant blue padded top
- Rubberised feet to prevent slipping
- Can be stacked with additional steps to increase height in 50mm (2") increments



Dimensions

- Height: 100mm (4")
- Width: 500mm (19½")
- Depth: 400mm (16")

TECHNICIAN

CEAx

Maximum User Weight: 127kg (20st)

N65731 - Portable Shower Screen

Description

- 3-panel portable, folding, shower screen
- For use in level access shower area

Dimensions

- Height: 815mm (32")
- Length when extended: 1474mm (58")
- Length when folded: 530mm (21")
- Product weight: 9kg (20lbs)



TECHNICIAN

N68185 - Savannah 2" Toilet Seat

Description

- Toilet seat for those requiring additional support and assistance sitting and standing.
- Gentle slope assists with lowering and raising.
- Contoured seat for comfort.
- Front and rear cut-outs aid personal cleansing.
- Fits to the toilet bowl with two non slip rubber buffers on L-shaped moulded plastic brackets, and a front lip

Dimensions

- Height: 50mm (2")
- Width: 370mm (14½")
- Length: 407mm (16")
- Aperture: 266x207mm (10½x8")



TECHNICIAN

CEAx

Maximum User Weight: 190kg (30st)

N68197 - Savannah 4" Toilet Seat

Description

- Toilet seat for those requiring additional support and assistance sitting and standing.
- Gentle slope assists with lowering and raising.
- Contoured seat for comfort.
- Front and rear cut-outs aid personal cleansing.
- Fits to the toilet bowl with two non slip rubber buffers on L-shaped moulded plastic brackets, and a front lip
- Removable.
- It fits most standard sizes of toilet.

Dimensions

- Height: 100mm (4")
- Width: 370mm (14½")
- Length: 407mm (16")
- Aperture: 266 x 207mm (10½ x 8")



TECHNICIAN

CEAx

Maximum User Weight: 190kg (30st)

N05266 - 6" Toilet Seat

Description

- 6" raised toilet seat
- Easy to fit and remove
- Secure by large fluted hand wheels
- Gentle slope to assist raising and lowering
- Rear cut out for comfort
- Front cut out to aid personal cleaning

Dimensions

- Height: 150mm (6")
- Width: 350mm (13 ¾")
- Length: 400mm (15 ¾")
- Aperture: 265mm x 210mm (10½" x 8¼")



TECHNICIAN

CEAx

Maximum User Weight: 190kg (30st)

M00869 - Toilet Frame

Description

- Height adjustable toilet frame
- Moulded arm supports for comfort
- Slip resistant rubber ferrules for extra safety

Dimensions

- Height: 635mm–785mm (25"-31")
- Width between arms: 430mm (17")
- Footprint: 502mm x 450mm (20" x 18")



Please state Height when prescribing this item.

TECHNICIAN

CEAx



Maximum User Weight: 190kg (30st)

M00870 - Toilet Frame Adjustable

Description

- Freestanding toilet frame
- Height and width adjustable
- Simple e-clip adjustment

Dimensions

- Height: 635mm–785mm (25–31")
- Width between arms: 420mm-550mm (16½"-21½")
- Footprint: 590mm–720mm x 470mm (23½"-28½" x 18½")



Please state Height when prescribing this item.

TECHNICIAN

CEAx



Maximum User Weight: 190kg (30st)

M24308 - Toilet Frame Fixing Kit

Description

- Toilet Frame Floor Fixing Kit to provide extra stability
- Requires a total of 12 screws to fix the frame to the floor (using fixing methods appropriate to the floor's construction)
- Use with M00869 - Toilet Frame & M00870 - Toilet Frame Adjustable



TECHNICIAN

CEAx

Maximum User Weight: 190kg (30st)

M66625 - Toilet Frame with Raised Seat

Description

- Freestanding, height adjustable toilet frame with moulded seat
- Seat has a gentle slope from back to front
- Front cut-out to aid personal cleaning
- Rubber ferrules to feet

Dimensions

- Height to Top of Seat: 400mm-550mm (15¾"-21½")
- Seat Width: 406mm (16")
- Seat Depth: 450mm (17¾")
- Height to Top of Frame: 610mm-770mm (24"-30")
- Aperture: 260 x 220mm (10" x 8½")
- Footprint: 520mm x 515mm-565mm (20½" x 20"-22")
- Width between arms: 450mm (17¾")



Please state Height when prescribing this item.

TECHNICIAN

CEAx



Maximum User Weight: 190kg (30st)

F18864 - Commode Height Adjustable

Description

- Adjustable height commode.
- Padded seat cover.
- Lift up toilet seat to fit and remove potty

Dimensions

- Seat height: 470mm–22mm (18½"-24")
- Width between arms: 490mm (19")
- Seat to top of back: 330mm (13")
- Seat to top of arms: 191mm (7½")



Please state Height when prescribing this item.

TECHNICIAN

CEAx



Maximum User Weight: 127kg (20st)

H87668 - Trolley

Description

- Enables the user to transport items around the home
- Not to be used as a mobility aid
- Height adjustable
- Lower shelf is set forward for ease of walking or dining
- High lipped trays prevent items from falling on the floor

Dimensions

- Height: 840mm-990mm (33"-39")
- Width: 460mm (18")
- Depth: 510mm (20")
- Weight: 7.5kg (17lb)
- Tray width: 410mm (16")
- Lower tray height: 330mm-480mm (13"-19")
- Upper tray height: 650mm-810mm (26"-32")
- Castors: 100mm (4")



Please state Height when prescribing this item.

TECHNICIAN

CEAx



Maximum User Weight: 177kg (28st)

L29877 - Buckingham Caddy

Description

- For use when someone is dependent on a walking frame that has been prescribed by a physiotherapist.
- To assist users of walking frames to independently move items around the home.
- Recommended for walking frames with 2 or 4 wheels.
- Supplied with a tray and mug holder.
- The person must be made aware they should not overload the caddy.
- One large and one small compartment with a tray to safely carry a plate and a mug holder plus anti-spill mug cover.

Dimensions

- Height: 200mm (8")
- Depth: 310mm (12¼")
- Width (with tray): 457mm (18")
- Weight: 1kg.
- Maximum carry weight: 2kg

**** The Buckingham Caddy requires a completed risk assessment / check list ****

**** TO BE USED WITH WHEELED WALKERS ONLY. ****

TECHNICIAN

CEAx



Maximum Carry Weight: 2kg

H87679 - Perching Stool

Description

- Enables people who are unable to stand for long periods of time to work in the kitchen, or wash in the shower and at the bathroom basin.
- Angled seat (front lower than the back) allows the User to rest then stand with minimal effort
- Height adjustable
- Splayed legs for stability
- Padded seat (foam covered with heat sealed PVC)
- Waterproof and wipe clean seat

Dimensions

- Height: 530mm-690mm (21"-27")
- Overall Width: 580mm (23")
- Seat Width: 340mm (13")
- Depth: 410mm (16")
- Seat Depth: 280mm (11")
- Weight: 4kg (8lb)



Please state Height when prescribing this item.

TECHNICIAN

CEAx



Maximum User Weight: 165kg (26st)

F19893 - Perching Stool with Arms & Back

Description

- Enables people who are unable to stand for long periods of time to work in the kitchen, or wash in the shower and at the bathroom basin.
- Angled seat (front lower than the back) allows the User to rest then stand with minimal effort
- Height adjustable
- Tubular arm rests
- Tubular back rest
- Soft padded vinyl seat
- Made from epoxy coated steel for strength and durability

Dimensions

- Seat Height: 515mm-640mm (20¼"-25¼")
- Seat Width: 345mm (13½")
- Seat Depth: 300mm (12")
- Footprint lowest: 510mm-370mm (20"x14½")
- Footprint highest: 550mm-400mm (21½"-15¾")



Please state Height when prescribing this item.

TECHNICIAN

CEAx



Maximum User Weight: 160kg (25st)

L33236 - Perching Stool with Arms & Padded Back

Description

- Enables people who are unable to stand for long periods of time to work in the kitchen, or wash in the shower and at the bathroom basin.
- Angled seat (front lower than the back) allows the User to rest then stand with minimal effort
- Height adjustable
- Tubular arm rests
- Padded back rest
- Soft padded vinyl seat
- Made from epoxy coated steel for strength and durability

Dimensions

- Seat Height: 515mm-640mm (20¼"-25¼")
- Seat Width: 345mm (13½")
- Seat Depth: 300mm (12")
- Padded Back: 300mm x 110mm x 35mm (11¾" x 4¼" x 1½")
- Footprint Lowest: 510mm x 370mm (20" x 14½")
- Footprint Highest: 550mm x 400mm (21½" x 15¾")



Please state Height when prescribing this item.

TECHNICIAN

CEAx

Clinical / Practical
Reasoning Required



Maximum User Weight: 160kg (25st)

N70854 - High Density 4" Foam Cushion

Description

- A high-density cushion designed to provide additional height to an existing chair or sofa that cannot be raised
- Ideally the cushion should be used beneath an existing chair cushion

Dimensions

- 460mm x 460mm x 100mm (18" x 18" x 4")



**** Please Note: This cushion has no pressure relieving properties ****

TECHNICIAN

CEAx

F19571 - Leg Lifter with Loop

Description

- Stiffened blue nylon strap
- Enables a User to move a lower limb when in a seated position
- Upper loop fits around the hand and wrist
- Lower loop fits around the foot

Dimensions

- Length (the rigid section): 1060mm (41½")
- Width of foot ring (Max): 130mm (5")
- Product weight: 87g (3oz)



TECHNICIAN

CEAx

M48192 - Bed Stick (Slatted or Divan)

Description

- For those who require assistance getting in and out of bed or for repositioning themselves once in bed
- Reversible plastic base allows use on both divan AND slatted beds
- Fits between mattress and bed base
- The slimmer base minimizes the impact of the product under the mattress
- Supplied with two easy to fit straps to hold securely in place

Dimensions

- Height from bed (slatted): 515mm (20")
- Width: 335mm (13")
- Length: 615mm (24")
- Base depth: 15mm (½")

**** Bed levers require a completed risk assessment and fitting guide ****

**** Not suitable for sprung based beds or profiling beds ****

TECHNICIAN

CEAx



Maximum User Weight: 160kg (25st)

M54234 - Bed Stick Double Ended (Divan)

Description

- For use on divan beds only
- For those who require assistance getting in and out of bed or for repositioning themselves once in bed
- Fits between mattress and bed base
- Width adjustable up to 1600mm beds
- One rail side can be reversed if rail only required on one side

Dimensions

- Width: 960mm-1600mm (38"-63")
- Depth: 356mm (14")
- Handle Height: 190mm (7½")
- Width of Loop: 356mm (14")
- Height of Loop: 432mm (17")

**** Bed levers require a completed risk assessment and fitting guide ****

**** Not suitable for profiling beds. Should only be considered for sprung based beds following a risk assessment due to the risk of vertical deflection causing instability to the service user when transferring ****



TECHNICIAN

Maximum User Weight: 190kg (30st)