A6 STANDING TRANSFER STAND & TURN FRAME - WITH COUNTERBALANCE

- 1
- Carer to ensure correct positioning of frame, and chair.
- Person to shuffle bottom to front of seat.
- Person to position feet on base plate.
- Carer to extend counterbalance bar underneath chair.
- Person to hold frame, lean forward with head raised and stand.
- 3
- Once person is standing, carer to retract counterbalance bar and swivel frame.

- 4
- Carer to extend counterbalance bar underneath chair.
- Person to sit down.
- Once seated, carer to remove frame.
- N.B. Once standing with yellow bar extended, carer no longer needs to counterbalance frame so may assist with personal care tasks.





